


# How to use face masks properly




## Wash hands first

Always make sure you wash your hands thoroughly before and after touching a mask.



## Masks

There are many different kinds of masks, depending on the task.



## The right side

There is a metal clip at the top of the mask.



## Placement

Place the metal clip across the top of your nose.



## Attach the mask

Attach the mask by pulling the elastic bands over your ears.



## Stretch down

Stretch the mask down, so that it covers your chin.



## Adjust

Bend the metal clip around your nose so that it sits securely.



## Taking off the mask

Pull the elastic bands away from your ears.



## Disposal

Always place the used mask in a closed rubbish bin for secure disposal.